


Universal Set Large II

With extra large drape for anaesthesia and extra large side drapes


Ref: 102-050-S

Component	Size (cm)	Material	Pcs
Table cover	150 x 200	2SBL	1
MAYO cover	75 x 145	PE70+NV28	1
Absorbent towel	32 x 32	MB40	4
Drape with adhesive edge	170 x 180	DS3, AP5	1
Drape with adhesive edge	180 x 300	DS2, AP5	1
Drape with adhesive edge	100 x 100	DS2, AP5	2
OP-tape	10 x 50	NV28, AP10	1
DINAFIX	3 x 20		1
Wrapping sheet	75 x 80	PF 2	1



ORL Drape


Ref	106-102-S
Size (cm)	200 x 340
Material	2 SBL
Aperture	∅ 7 cm with adhesive
Reinforced area	60 x 75 made from NV28



Ear Surgery Set

Ref: 106-002-S


Component	Size (cm)	Material	Pcs
Wrapping sheet	100 x 100	PF2	1
Table cover	100 x 150	2 SBL	1
MAYO cover	75 x 145	PE, NV28	1
Surgical drape Waterproof	150 x 200	2 SBL, AP5	1
Self-adhesive pouch small	22 x 28	PE-foil, AP2,5, FX	1
Surgical drape	150 x 155	SMS-blue	2
Surgical drape	51,6 x 80	SMS-blue	1
Protective gown Basic	XL		2
Self-adhesive pouch large	28 x 44	PE-foil, FX, AP5, omnifilm	1
Protective gown Basic	L		1
Fenestrated drape	90 x 90	DS2, AP15	1



Strumectomy Set


Ref: 106-003-S

Component	Size (cm)	Material	Pcs
Wrapping sheet	100 x 100	PF2	1
Table cover	100 x 150	2 SBL	1
MAYO cover	75 x 145	PE, NV28	1
Surgical drape	150 x 155	SMS-blue	2
Protective gown Basic	XL		2
Surgical gown Classic	L		1
Protective gown Basic	L		2
Self-adhesive pouch small	22 x 28	PE-foil, FX, AP2,5	1
Self-adhesive pouch large	28 x 44	PE-foil, FX, AP5, omnifilm	2
Struma drape	150 x 300	2 SBL, DS2, IF	1



Tracheostomy Drape

Ref	501-008-S
Size (cm)	155 x 300
Material	SMS
Aperture	∅ 10 cm with adhesive
Reinforced area	75 x 100 made from DS2



Ophthalmic Drape

With integrated self-adhesive incision foil

△ Fluid collection pouches with plastic brace at both sides of IF aperture

Ref	201-101-S
Size (cm)	155 x 122
Material	SMS, 2x60-0101, FX
Aperture	Pentagonal 11 x 12 cm with IF

