


Hip TEP Set

Ref: 301-450-S


Component	Size (cm)	Material	Pcs
Wrapping sheet	75 x 80	PF2	1
Table cover	150 x 200	2 SBL	1
Drape 2 SBL	150 x 200	2 SBL	1
Drape DS2 with adhesive edge	150 x 240	DS2, AP5-80	1
Hip U-drape combined	200 x 280	DS3, DS2 - laminate, AP15	1
OP-tape	10 x 50	NV28, AP10	2
Limb cover - large	40 x 120	DS2	1
MAYO cover reinforced	75 x 145	PE, DS2	1



Head Surgery Set

Ref: 402-106-S


Component	Size (cm)	Material	Pcs
Table cover	150 x 200	2SBL	1
MAYO cover reinforced	75 x 145	PE70 + DS2	1
Absorbent towel	32 x 32	MB40	4
Head surgery drape	220 x 240	DS3	1
OP-tape	10 x 50	NV28, AP10	1
DINAFIX	3 x 20		1
Wrapping sheet	75 x 80	PF2	1



Spine Surgery Set


Ref: 401-116-S

Component	Size (cm)	Material	Pcs
Table cover	150 x 200	2SBL	1
MAYO cover reinforced	75 x 145	PE70 + DS2	1
Absorbent towel	32 x 32	MB40	4
Spine surgery drape	150 x 280	DS3, IF	1
F/C pouch	28 x 44	PE50	1
OP-tape	10 x 50	NV28, AP10	1
DINAFIX	3 x 20		1
Wrapping sheet	75 x 80	PF2	1




Head Surgery Drape

Ref	401-201-S
Size (cm)	220 x 240
Material	DS3
Aperture	10 x 20 with IF
F/C pouch	28 x 44



Spine Surgery Drape

Ref	401-204-S
Size (cm)	255 x 315
Material	DS3
Aperture	9 x 20 cm with adhesive



Angiography Drape

With double transparent sides and highly absorbent reinforced area. Equipped with two adhesive femoral oval fenestration with incision foil

Ref	501-115-S
Size (cm)	220 x 365
Material	DS2, PE 70 µm
Aperture	∅ 12 cm / drape, ∅ 8 cm with IF
Reinforcement	65 x 76 cm / DS2 Lam

