


TUR Drape combined


Ref	603-205-S
Size (cm)	220 x 300
Material	DS3, PE blue
Aperture	20 x 20 with filter



TUR Set

Ref: 603-107-S

Component	Size (cm)	Material	Pcs
Table cover	100 x 150	2SBL	1
Assistant protective gown Basic	L		1
Urology gown	XXL		1
TUR drape	150 x 220	DS2 + SMS	1
Limb cover	78 x 100	SMS	2
Wrapping sheet	75 x 80	PF2	1



Gyneacology Drape


Ref	603-206-S
Size (cm)	75 x 175
Material	DS2
Aperture	10 x 14



Perineal Set

Ref: 604-001-S


Component	Size (cm)	Material	Pcs
Wrapping sheet	100 x 100	PF2	1
Table cover	100 x 150	2 SBL	1
Drape DINASAFE with adhesive edge	50 x 77,5	SMS-blue, AP2,5	1
Laparotomy drape III short	180 x 200	2 SBL, NV28, IF45x50, AP5, AP2,5	1
Perineal drape III	75 x 120	DS2, AP15, AP2,5	1
OP-tape	10 x 50	NV28, AP10	3
Limb cover III	77,5 x 120	SMS-blue	2
Armamentarium drape	50 x 51,6	SMS-blue	1
Table cover	100 x 150	2 SBL	2
MAYO cover reinforced	75 x 145	PE, DS2	1
Surgical gown Classic	L		1
Absorbent towel	32 x 32	MB40	4



Rectum Surgery Set

Ref: 604-003-S


Component	Size (cm)	Material	Pcs
Wrapping sheet	75 x 80	PF2	1
Table cover	100 x 150	2 SBL	1
MAYO cover reinforced	75 x 145	PE, DS2	1
Limb cover	77,5 x 100	SMS-blue	2
TUR drape	150 x 260	SMS-blue, DS2	1



Gynaecology Set

Ref: 601-103-S

Component	Size (cm)	Material	Pcs
Wrapping sheet	75 x 80	PF2	1
Table cover	100 x 150	2 SBL	1
Combined drape with adhesive edge II	180 x 200	2 SBL, NV28, AP5, AP2,5	1
Perineal drape III	75 x 120	DS2, AP15, AP2,5	1
Limb cover III	77,5 x 120	SMS-blue	2
Table cover	100 x 150	2 SBL	1
Armamentarium drape	50 x 51,6	SMS-blue	1
OP-tape	10 x 50	NV28, AP10	2
MAYO cover	75 x 145	PE, NV28	1
Surgical gown Classic	L		1
Absorbent towel	32 x 32	MB40	4



Gynaecology Set for Laparoscopy

Ref: 601-132-S

Component	Size (cm)	Material	Pcs
Wrapping sheet	75 x 80	PF2	1
Table cover	150 x 200	2 SBL	1
Absorbent towel	32 x 32	MB40	2
Surgical drape	75 x 90	DS2	1
Limb cover	77,5 x 120	SMS	2
Patient drape	200 x 180	2 SBL., NV28, AP15, AP2,5	1
MAYO cover	75 x 145	PE, DS2	1

