


Surgical Gown OPTIMUM

- △ Made from non-woven textile hydrophobic material SMS that protects against fluid penetration
- △ Gowns are fixed with four ties system
- △ Sleeves are finished with knitted elastic cuffs
- △ Surgical gowns are a tie-round type, which enables to keep the sterility of surgeons back side
- △ **Fulfil requirements of ČSN 13795 standard**

Size	Lenght in cm	Girth in cm	Ref
M	125	150	353023-S
L	135	155	353033-S
XL	145	164	353043-S
XXL (extra long)	155	164	353053-S




Surgical Gown OPTIMUM

With reinforcement

- △ Gown is reinforced at front side and sleeves with impermeable material that assures perfect protection of surgeon during major surgeries
- △ **Fulfil requirements of ČSN 13795 standard**


Size	Lenght in cm	Girth in cm	Ref
M	125	150	353123-S
L	135	155	353133-S
XL	145	164	353143-S
XXL (extra long)	155	164	353153-S



Surgical Gown CLASSIC

- △ Made from non-woven textile hydrophobic material SMS that protects against fluid penetration
- △ Gowns are fixed with four ties system
- △ Sleeves are finished with knitted elastic cuffs
- △ Surgical gowns are a tie-round type, which enables to keep the sterility of surgeons back side
- △ Antistatic and alcohol repelent modification
- △ **Fulfil requirements of ČSN 13795 standard**

Size	Lenght in cm	Girth in cm	Ref
M	125	150	351023-S
L	135	155	351033-S
XL	145	164	351043-S
XXL (extra long)	155	164	351053-S




Surgical Gown CLASSIC

With reinforcement

- △ Gown is reinforced at front side and sleeves with impermeable material that assures perfect protection of surgeon during major surgeries
- △ **Fulfil requirements of ČSN 13795 standard**


Size	Lenght in cm	Girth in cm	Ref
M	125	150	351123-S
L	135	155	351133-S
XL	145	164	351143-S
XXL (extra long)	155	164	351153-S



Surgical Gown PREMIUM

- △ Made from hydrophobic non-woven material Sontara with enhanced air permeability. Protects from fluid penetration
- △ Gowns are fixed with four ties system
- △ Sleeves are finished with knitted elastic cuffs
- △ Surgical gowns are a tie-round type, which enables to keep the sterility of surgeons back side
- △ Recommended for long time procedures
- △ Very comfortable to wear
- △ **Fulfil requirements of ČSN 13795 standard**

Size	Lenght in cm	Girth in cm	Ref
M	125	150	350023-S
L	135	155	350033-S
XL	145	164	350043-S
XXL (extra long)	155	164	350053-S



Surgical Gown PREMIUM

With reinforcement

- △ Gown is reinforced at front side and sleeves with impermeable material that assures perfect protection of surgeon during major surgeries
- △ **Fulfil requirements of ČSN 13795 standard**

Size	Lenght in cm	Girth in cm	Ref
M	125	150	350123-S
L	135	155	350133-S
XL	145	164	350143-S
XXL (extra long)	155	164	350153-S

